

HPV Vaccine Plan

INTRODUCTION

In VIP, the new HPV (human papilloma virus) Vaccine Plan will be implemented similar to MeNZB plan.

In the entire plan, there are 2 types of vaccines to be given :


1. Human Papilloma virus vaccine, bivalent = HPV bivalent
2. Human Papilloma virus vaccine, quadrivalent = HPV quadrivalent

Start from 1 September 2008, as per request by Ministry - the plan is implemented only for HPV quadrivalent. It is restricted to females born on or after 01/01/1990.

HPV Plan Installation

To include the HPV Plan into the Immunisation Plan setup, you need to install the plan.

Go to File > Utility > HPV Vaccination Programme.



Choose the option **HPV quadrivalent vaccine only**.

Click 'Activate' button to install the plan. Normally it will take not more than a minute.

The installation needs to be done only ONCE.

'Deactivate' button is provided in case user install the plan by mistake before the time.

If there is patient has been given the vaccines in the plan in the immunisation records, it will put the Termination date in.

If there is NO patient has been given, it will delete the plan in the Setup.

Fee Setup

Once HPV plan has been installed, the fee needs to be set up.


Go to **Setup > Medical > Immunisation**.

Highlight **HPV quadrivalent** from the list, then press **ENTER** or click **MODIFY** button.

Go to **Fees** tab.

Fill in the **Subsidy** field with the current subsidy fee, \$ 20.16.

Then press **TAB**. It will update **Total Fee** field automatically.



The screenshot shows a software window with three tabs: 'Browse', 'Details', and 'Fees'. The 'Fees' tab is active. Inside the window, there are four input fields: 'Subsidiser' (a dropdown menu showing 'IMM Claiming'), 'Subsidy' (a text box with '\$20.16'), 'Surcharge' (a text box with '\$0.00'), and 'Total Fee' (a text box with '\$20.16'). Below these fields is a button labeled 'Review Charges Setup'. At the bottom of the window, there are four navigation arrows (back, left, right, forward), a 'New Record' button, and 'OK' and 'Cancel' buttons.


Press **F12** or click **OK** button to save the changes.

This setup needs only to be done ONCE.

Immunisation Screen (Shift+F9)

On the below part of the Shift+F9 screen, there is a tick box for HPV Plan which works exactly similar to the MeNZB plan.

Once it's ticked, all of the records of this plan will be messaged to the NIR regardless the patient casual status and the opt off status.



The screenshot shows a section titled 'Plans Enrolled for:'. It contains seven checkboxes with labels: 'Standard Plan' (checked), 'MeNZB Plan', 'BCG Plan', 'HepB Positive Mother Plan', 'Pneumococcal Plan', 'Pre/Post Splenectomy', and 'HPV Plan'. The 'HPV Plan' checkbox is circled in red.

You can remove a scheduled HPV plan (provided that the immunisation has not been given) by clicking on the HPV Plan box to remove the tick.

Please note that HPV Plan tick box is only enable for Female patients.

And the schedule of the HPV Plan is determined by the age at dose one.
 Details of the plan are given in the table below.

Age of child	Vaccine	Dose number	Indicator
Age dose 1 due	HPV quadrivalent	1	5 = Primary course
Age dose 1 due + 2 months	HPV quadrivalent	2	5 = Primary course
Age dose 1 due + 6 months	HPV quadrivalent	3	5 = Primary course

For the first time once you tick the HPV Plan, all the vaccines in the HPV plan with 'Unknown' in the Age Due.

Age Due	Vaccination	Status	Date Due	Date Done
Unknown	HPV, quadrivalent 1		//	
Unknown	HPV, quadrivalent 2		//	
Unknown	HPV, quadrivalent 3		//	

Once you administer the first HPV, bivalent - the program will calculate the due dates of the remaining vaccines. Please see the below screen.

Age Due	Vaccination	Status	Date Due	Date Done	Batch	Site
18y+2m	HPV, quadrivalent 1	✓ DONE at 18y+2m	20/08/2008	20/08/2008	BN8976	Deltoid, I/Mu
18y+4m	HPV, quadrivalent 2	Due in 61 days	20/10/2008			
18y+8m	HPV, quadrivalent 3		18/02/2009			

For patients opting off please refer to the process given to you by the NIR or contact your DHB Administrator.

Our Support team have created a report with the following criteria to enable you to import and run for a mail merge:
 The report is called HPV - and will give you a list of all patients born between 01/01/90 to